

Galaxy 2000 Alternating Air Mattress



A robust, reliable mattress suitable for home, nursing home and hospital use.

The Galaxy 2000 is a full replacement therapy mattress offering proven, effective, pressure ulcer treatment. It also offers prevention for users who are at a very high risk level or, have existing open tissue damage.
EPUAP Ulcer Grade I to IV

Features

Offers a 60 degree incline with minimal or no cell compression, allowing patients to sit up comfortably with minimal shear.

A foam under layer to provides stability when manoeuvring or getting in and out of bed.

Adjustable to patient's weight for optimum pressure relief.

Waterproof, vapour permeable, two way stretch covers for optimum patient comfort and care.

Sealed top surface seams and full depth side flaps minimise the risk of liquid ingress.

Power unit is both energy efficient and quiet in operation.

Can be set to stay inflated for approximately 24 hours incase of mains power failure.

Made to measure sizes available including Bariatric sizes.

Clinical Factors

Risk Level: Very high risk

Application: Up to and including category IV

Weight limit: 280kg/45 stone

Inflated cell height: 115mm

Cycle time: 15 minutes

