

Stellan Complex Care System

Comfort and safety for those living with a range of care issues including Huntington's and Dementia



HUNTINGTON'S

DEMENTIA

PARKINSON'S

ALZHEIMER'S

- The Stellan System is made with strong padded sides equipped to absorb the force of involuntary movements.

- Sides can be made up to 700mm high from the top of the mattress which will reduce the risk of falls.

- Optional split sides are available. These allow the person in the bed to be attended to and still be safe from rolling out of the bed. This provides effective protection against injury to patient and carer.

- Electric height adjustment provides safe and comfortable access when transferring in and out of the bed.

- A base sheet attached to the four padded sides helps reduce entrapment issues.


- Profiling the head section may help respiration. Applying the knee brake can help to prevent the person slipping down the bed when experiencing reflexive twisting and arching movements.

- For someone with full mobility at risk of falling/climbing out of bed a cut out side can be a safe option. The part side will provide protection from rolling out but won't restrict the person from exiting the bed feet first. This is also available on an extra low bed.

- There are many options and because everyone's issues, preferences and needs are different we always offer a home visit to assess the most suitable option from one of our product specialists.



TRENDELENBURG TILT



SPLIT SIDE



CUT OUT SIDE



Kinderkey Healthcare Ltd. Unit 1, Cilmedw Way, Llangollen, LL20 8AG

Tel: 01978 820714 contact@kinderkey.co.uk www.kinderkey.co.uk